EL

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REF 871
Size: 22in x 17.25in (55.9cm x 43.8cm)

MAXI-THERM® LITE LOCALIZED THERAPY PAD Thermal Regulating System

WARNINGS

- Licensed healthcare practitioner’s order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a licensed healthcare practitioner. Higher risk patients; including pediatric, temperature sensitive and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to bony prominenc-es. Notify licensed healthcare practitioner for changes in patient condition.
- Follow licensed healthcare practitioner order/hospital policy for discontinuation of therapy. Notify licensed healthcare practitioner for changes in patient condition.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Leaky pads or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not place additional heat sources between the patient and Maxi-Therm® Lite. Tissue damage may result.
- Do not block the fluid pathways with positioning accessories. Patient therapy may be interrupt-ed.
- Do not place several layers of sheets or other materials that will reduce the heat transfer between the blanket and the patient.
- Do not place sharp objects in contact with the Maxi-Therm® Lite pad as punctures may occur. Do not use pins to secure the pad.
- Do not use the device distal to arterial cross clamping. Non-observance can lead to thermal injury.
- The Maxi-Therm® Lite pad is single patient use only. Reuse may result in biocontamination or cross-infection.

PRECAUTIONS

- Caution: Federal law restricts the device to sale by or on the order of a licensed healthcare practitioner.
- Presence of water or gel on the pad while the pad is not operating may cool the patient.
- Use of unapproved equipment or hoses is not recommended. The Maxi-Therm® Lite has been designed and tested for use with Gentherm thermal regulating systems only.
- If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury
- Do not place bean bags or other firm positioning devices under the blankets/pads. Skin damage may occur.
- In surgery, heating may affect toxicity of certain prep solutions which may cause skin injury when remaining on patients’ skin during application of prolonged heat.

CE 0344

GENTHERM

Gentherm Medical, LLC
12011 Mosteller Road
Cincinnati, Ohio, USA 45241
Tel: 1-800-989-7373
www.gentherm.com

- Keep the area between the patient and the blanket dry. Excessive moisture may result in skin damage.

INDICATIONS FOR USE

The Maxi-Therm® Lite hyper-hypothermia pad is intended to be used by trained healthcare professionals as part of Gentherm’s thermal regulating systems, specifically the Blanketrol® II, Blanketrol® III, Micro-Temp® LT, or Norm-O-Temp®, to lower or to raise a patient’s temperature and/or maintain a desired patient temperature through conductive heat transfer.

CONTRAINDICATIONS