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
Size: 24in x 60in (61cm x 152.4cm)

WARNINGS

- Licensed healthcare practitioner’s order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a licensed healthcare practitioner. Higher risk patients; including pediatric, temperature sensitive patients and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to bony prominenc-es. Notify licensed healthcare practitioner for changes in patient condition.
- Follow licensed healthcare practitioner order/hospital policy for discontinuation of therapy. Notify licensed healthcare practitioner for changes in patient condition.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Leaky blankets or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not place additional heat sources between the patient and Maxi-Therm®. Tissue damage may result.
- Do not block the fluid pathways with positioning accessories. Patient therapy may be interrupt-ed.
- Do not place several layers of sheets or other materials that will reduce the heat transfer between the blanket and the patient.
- Leaky blankets or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not place sharp objects in contact with the Maxi-Therm® blanket as punctures may occur. Do not use pins to secure Maxi-Therm® blanket.
- Do not use the device distal to arterial cross clamping. Non-observation can lead to thermal injury.
- The Maxi-Therm® blanket is intended for single patient use only. Reuse may result in biocontamination or cross-infection.

- Keep the area between the patient and the blanket dry. Excessive moisture may result in skin damage.

INDICATIONS FOR USE

The Maxi-Therm® hyper-hypothermia blanket is intended to be used by trained healthcare professionals as part of Gentherm’s thermal regulating systems, specifically the Blanketrol® II, Blanketrol® III, Hemotherm®, or Norm-O-Temp®, to lower or to raise a patient’s temperature and/or maintain a desired patient temperature through conductive heat transfer.

CONTRAINDICATIONS

None known.

SET UP

- Unfold blanket on bed, stretcher, or operating room (OR) table with hoses routed without kinks towards the unit. Check blanket before each use for product integrity.
- Check the equipment reservoir to

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