ΕΛ

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**REF 193**  
**Size:** 12in x 18in (30.5cm x 45.7cm)

**WARNINGS**

- Licensed healthcare practitioner’s order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a licensed healthcare practitioner. Higher risk patients; including pediatric, temperature sensitive and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to bony prominences. Notify licensed healthcare practitioner for changes in patient condition.
- Follow licensed healthcare practitioner order/hospital policy for discontinuation of therapy. Notify licensed healthcare practitioner for changes in patient condition.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Leaky blankets or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not place additional heat sources between the patient and PlastiPad®. Tissue damage may result.
- Do not block the fluid pathways with positioning accessories. Patient therapy may be interrupted.
- Do not place several layers of sheets or other materials that will reduce the heat transfer between the blanket and the patient.
- Do not place sharp objects in contact with the PlastiPad® blanket as punctures may occur. Do not use pins to secure PlastiPad® blanket.
- Do not use the device distal to arterial cross clamping. Non-observation can lead to thermal injury.

**PRECAUTIONS**

- Caution: Federal law restricts the device to sale by or on the order of a licensed healthcare practitioner.
- Presence of water or gel on the blanket while the blanket is not operating may cool the patient.
- Use of unapproved equipment or hoses is not recommended. The PlastiPad® has been designed and tested for use with Gentherm thermal regulating systems only.
- If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury.
- Do not place bean bags or other firm positioning devices under the blankets/pads. Skin damage may occur.
- In surgery, heating may affect toxicity of certain prep solutions which may cause skin injury when remaining on patients’ skin during application of prolonged heat.
- Keep the area between the patient and the blanket dry. Excessive moisture may result in skin damage.

**EC REP**  
**CE0344**  
Authorized European Representative  
CEpartner4U, BV  
Esdsmoerlaan 13  
3951 DE Maam  
The Netherlands  
www.CEpartner4U.com

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